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
期刊文献 Certification of Beryllium Mass Fraction in SRM 1877 Beryllium Oxide Powder Using High-Performance Inductively Coupled Plasma Optical Emission Spectrometry with Exact Matching

作者: [Winchester, Michael R. mrw@nist.gov](mailto:mrw@nist.gov) [Turk, Gregory C.](#) [Butler, Therese A.](#) [Oatts, Thomas J.](#) [Coleman, Charles](#) [Nadratowski, Donald](#) [Sud, Ritu](#) [Hoover, Mark D.](#) [Stefaniak, Aleksandr B.](#)

High-performance inductively coupled plasma optical emission spectrometry (HP-ICP-OES) was used to certify the Be mass fraction in National Institute of Standards and Technology (NIST) Standard Reference Material (SRM) 1877 Beryllium Oxide Powder. The certified value and expanded uncertainty expressed at a 95% confidence level is **(0.3576 ± 0.0024) g/g.**

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- **The Customer Request**

- Demonstrate that the Optima is capable of analyzing high accurate and precise determinations of PMGs

- **The Elements**

- Au, Pt, Pd, Rh, Ir, Ag, Ru
- e.g. **ISO 11494 & 11495**

- **The Problem**

- Analytical Precision Relatively Poor : 0.20 - 0.40%

$$\bar{m}_{Pt,fin} = \frac{1}{5} \left(\sum_{n=1}^5 m_{Pt} \right)$$

The RSD of m_{Pt} shall not exceed 0.30 %.

INTERNATIONAL
STANDARD

ISO
11494

Second edition
2014-12-01

Jewellery — Determination of
platinum in platinum jewellery alloys
— ICP-OES method using yttrium as
internal standard element

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- NIST produces and certifies 70 Spectrometric Solution SRMs
 - 1000's sold annually for calibration of ICP, AAS
 - typically 10 mg/g
 - **~0.3% uncertainty**
- transfer standards used for calibration

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Anal. Chem. 2001, 73, 4821–4829

Single-Element Solution Comparisons with a High-Performance Inductively Coupled Plasma Optical Emission Spectrometric Method

Marc L. Salit,* Gregory C. Turk, Abigail P. Lindstrom, Therese A. Butler, Charles M. Beck II, and Bruce Norman

Chemical Science and Technology Laboratory, National Institute of Standards and Technology, 100 Bureau Drive, Gaithersburg, Maryland 20899-8391

A solution-based inductively coupled plasma optical emission spectrometric (ICP-OES) method is described for elemental analysis with relative expanded uncertainties on the order of 0.1% relative. The single-element determinations of 64 different elements are presented, with aggregate performance results for the method and parameters for the determination of each element. The performance observed is superior to that previously reported for ICP-OES, resulting from a suite of technical strategies that exploit the strengths of contemporary spectrometers, address measurement and sample handling noise sources, and permit rugged operation with small uncertainty. Taken together, these strategies constitute high-performance ICP-OES.


quantitatively different from that previously reported for ICP-OES results. Different chemical information becomes available with this ability to discriminate between solutions of very similar composition—an improvement in concentration resolution. Additionally, at this level of uncertainty, uncertainty is dominated by effects other than spectroscopic measurement (e.g., sample handling and preparation).

Measurements are performed using unmodified commercially available equipment and are significantly less costly than classical analysis. The opportunity to achieve results of the quality expected of classical analysis with instrumental efficiency enables new analytical applications. In contrast to classical analysis, instrumental automation permits the analysis of multiple samples with little incremental cost per sample. Similar to classical analysis, this method is suitable for determination of “major” and “minor”

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Drift is major
source of
uncertainty

Drift Correction

$$S_{measured} = S_{truth} + \epsilon_{drift} + \epsilon_{noise}$$

- separate variables to estimate drift as $f(t)$
 - use grand mean for sample to estimate S_{truth}
 - create time series of deviations from mean
 - use **ALL** standards and samples
 - this is our estimate of $(\epsilon_{drift} + \epsilon_{noise})!$
 - use relative deviations for multiplicative drift, absolute deviations for additive drift
 - fit smooth function to these deviations
 - this is our estimate of $\epsilon_{drift}!$

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0584-8547/83 \$3.00 + .00
© 1983 Pergamon Press Ltd.

Improved performance using internal standardization in inductively-coupled plasma emission spectroscopy

S. A. MYERS and D. H. TRACY

Perkin-Elmer Corporation, Main Avenue, Norwalk, CT 06856, U.S.A.

(Received 30 June 1982, in revised form 25 February 1983)


Abstract—Significantly improved performance in ICP emission spectroscopy is possible using internal standardization given the proper choice of operating conditions. A dual monochromator detection system was used to study the correlation between two emission signals recorded simultaneously as several ICP parameters were varied. The behavior of the signal formed by dividing the analyte and reference signals was also measured.

The noise behavior of twenty analytically important elements was compared to that of a manganese ion line used as the internal standard for a single, compromise choice of ICP operating parameters (RF power = 1250 W, carrier gas flow rate = 0.7 or 0.9 l/min, and viewing height = 20 mm). The elements differed widely in excitation and ionization energies. Sample concentrations were chosen so that the analyte emission signals were sufficiently intense that the noise and drift contributions due to plasma background and shot noise could be ignored. Under these conditions very good signal correlation ($r \geq 0.95$) and similar noise proportionality factors were observed leading to noise reductions of $\geq 10 \times$ after signal division in almost every instance. After signal division the noise, defined as the relative standard deviation of the emission signal, (time constant = 1.0 s, $n = 20$) was less than 0.1%. The closely similar noise behavior is attributed to sample density fluctuations arising from corresponding fluctuations of the sample aerosol density. Good correlation is also observed between fluctuations in the plasma background simultaneously detected at two widely separated wavelengths when appropriate steps are taken to reduce shot noise.

At carrier gas flow rates above 1.0 l/min and at viewing heights below 15 mm or above 20 mm, the signal correlation is reduced or absent, yielding little or no improvement using internal standardization.

The use of internal standardization leads to other improvements in analytical performance in addition to noise reduction. Drift is reduced, and nebulizer performance is improved in that it is less sensitive to variations in liquid sample uptake rate, or nebulizer instabilities associated with high solids content sample introduction. The divided signal takes much less time to reach equilibrium following introduction of a new sample.


- precision enhancement
 - matched internal standard
 - can be ~20x improvement

Establishing Traceability to NIST SRMs
with High Performance ICP-OES
Marc Salit, NIST and Sue Evans, VHG


Results: Stability

- Monitor internal standards to check for instrument drift
 - Sc in radial mode, Y in axial mode

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

Make a simple Argon diagnostic method and monitor the RSD.

Run a 2v/v% nitric acid solution with four different conditions:

1. Pump on and nebulizer on
2. Pump on and nebulizer off
3. Pump off and nebulizer on
4. Pump off and nebulizer off

A	B	C	D	E	F	G	H	I
Exp_5 Ar-diagnose test								
Plasma power:	1350		Nebulizer type:	CF	resolution	normal		
Nebulizer flow:	0.85		Nebulizer tubing:	B/B	profiling	no		
Auxiliary flow:	0.2		Injector diameter:	2.0 mm	Replicates	10		
Plasma flow:	15		Peak Algorithm:	Area				
Pump speed:	1.5		Points per peak:	3	torch position	-3		
methode:	Exp_5		Results file:	Exp_5				
Read parameters	min 1	max 5						
test	Ar 420.069 % RSD							
pump on neb on	0.28							
pump on neb uit	0.18							
pump off neb on	0.14							
pump off neb out	0.03							

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Sequence No.: 1
Sample ID: Sample002
Analyst:
Logged In Analyst (Original) : lenovo
Initial Sample Wt:
Dilution: 1144X
Wash Time (before sample):

Autosampler Location:
Date Collected: 2017/10/18 16:47:10
Data Type: Reprocessed on 2017/10/19 14:03:11

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
Sequence No.: 3
Sample ID: Std2
Analyst:
Logged In Analyst (Original) : ICP1
Initial Sample Wt:
Dilution:
Wash Time:

Autosampler Location: 3
Date Collected: 12/18/2014 2:02:20 PM
Data Type: Reprocessed on 12/30/2014 3:08:56 PM

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

US — Swiss Zn Comparison

- comparison of NIST SRM against...
 - EMPA distilled Zn
 - high purity metals
- calibration curve for 20 solutions
 - 800 min expt
 - 0.0082% solution-to-solution variability

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Table 7 Precisions ($n = 10$) and recoveries of the method


Component	Found $w / (\%, \text{ mass fraction})$	Average $w / (\%, \text{ mass fraction})$	RSD/ %	Added $w / (\%, \text{ mass fraction})$	Found $w / (\%, \text{ mass fraction})$	Recovery/ %
Nd	29.90, 29.88, 29.84, 29.90, 29.87, 29.88, 29.85, 29.86, 29.81, 29.85	29.86	0.09	30.0	58.1	97
Fe	69.05, 69.07, 69.11, 69.05, 69.08, 69.07, 69.10, 69.09, 69.14, 69.09	69.08	0.04	70.0	142.95	103
B	1.051, 1.050, 1.045, 1.050, 1.051, 1.047, 1.053, 1.046, 1.053, 1.055	1.050	0.31	1.0	2.01	98

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



Stainless Steel Type 304L B.S. 81G Values expressed as weight%

Element and Wavelength	Reference Value	Avio500 Value
B 249.772	0.0005 ±0.001	0.00049 ±0.000005
Co 228.616	0.29 ±0.01	0.296 ±0.003
Cr 283.563	18.65 ±0.08	18.73 ±0.075
Cu 324.752	0.54 ±0.02	0.543 ±0.004
Mn 294.920	1.69 ±0.03	1.72 ±0.013
Mo 203.845	0.69 ±0.02	0.685 ±0.006

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
The screenshot shows a software window titled "Wavelength Table". It has a search interface with the following parameters:

- Center wavelength: 221.803 nm
- Wavelength range (+/-): 0.01 nm
- Wavelengths: All, Recommended
- Elements to include: Copper (Cu) (selected from a dropdown), All (button)
- Search button




Below the search parameters is a table with the following data:

	Elem.	State	Wavelength (nm)	Pref. Order	BEC (mg/L)	DL mg/L (W, P, &F)	Rel Sens(B)	Sig/Noise (S & T)	Intensity (W)
1	Cu	II	221.811		0.59	0.0170	32000	61.9	120.0
2									
3									

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

Sequence No.:	1	Autosampler:	
Sample ID:	200fe	Date Collecte	
Analyst:		Data Type:	Re
Logged In Analyst (Original):	chiss	Initial Samp	
Initial Sample Wt:		Sample Prep	
Dilution:			

Analyte	Mean Corrected Intensity	Calib. Conc. Units	Std.Dev.
Y 371.029	3182157.5	97.10 %	0.529
Co 238.892+	-3258283.5	-33.29 mg/L	0.255
Fe 238.204	27720743.8	175.9 mg/L	0.67
Sn 235.485+	568162.7	53.31 mg/L	0.007

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
Method Editor : UDA_Simple_X2

Define Elements

Method Description: UDA Simple Method

	Symbol	Wavelength (nm)	Name	Function
1	Co	238.892	Co 238.892	Analyte
2	Co	228.616	Co 228.616	Analyte
3	Fe	273.955	Fe 273.955	Analyte
4	Fe	238.204	Fe 238.204	Analyte
5	Sn	242.170	Sn 242.170	Analyte
6	Sn	235.485	Sn 235.485	Analyte
7	Y	371.029	Y 371.029	Int. Std.
8				
9				

Elements and wavelengths can be selected by clicking on one of the buttons to the right.....

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
Sequence No.: 1
Sample ID: 200fe
Analyst:
Logged In Analyst (Original) : chiss
Initial Sample Wt:
Dilution:

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The screenshot shows the 'Map Editor - For Element La' window. It contains a table with the following data:

Wavelength & Order	Wavelength Range	Column	Row	Peak Status	State	Relative Intensity	Available for Analysis	Auto Select for New Method
333.749 (101)	High	280.94	148.14	Not Located		1500000	<input checked="" type="checkbox"/>	<input type="checkbox"/>
379.478 (089)	High	375.15	113.55	Not Located		1200000	<input checked="" type="checkbox"/>	<input type="checkbox"/>
408.672 (082)	High	2.80	95.18	Not Located		1800000	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
412.323 (082)	High	429.22	96.09	Not Located		1500000	<input checked="" type="checkbox"/>	<input type="checkbox"/>
419.665 (080)	High	88.50	90.35	Not Located		500000	<input checked="" type="checkbox"/>	<input type="checkbox"/>

Below the table, there are tabs for 'Peak Location', 'Subarray', and 'General'. A message at the bottom says 'Indicates that we recommend you to use this new line.' There are 'Add Lines...' and 'Delete...' buttons.

An error dialog box is overlaid on the bottom right, with a red border around the message: 'Wavelength 398.852 is not available for this instrument. You are not able to add this line to the Map database.' The dialog has an 'OK' button.